

4 Select appropriate repositioning/refocusing options.

Select the channel used to guide the repositioning/refocusing process. A typical choice is the bottom bead image channel.

Refined recovery options

AcTrM - XYZ recovery options: X

Channel for rough recovery:

Phase 563-9 492-18

- Perform refined recovery

- Uses limited field-of-view for refined recovery

Channel for refined recovery:

Phase 563-9 492-18

- start the focus search by raising the objective

OK

Select to perform refined recovery

Select to limit refined recovery process to a limited field-of-view

Select to define upward motion as the first step of the focus search procedure