


5 *An "analysis" folder will appear in the file manager and contain a "p0" folder and an "analysisChoices.txt" file.

p0



De-drifted cell images

Top bead images from AcTrM


Contains flags for instances where analysis is skipped

De-drifted top bead images


Bottom bead images from AcTrM

Cell images from AcTrM

Amount of drift detected in each frame



analysisChoices.txt*



Channels

Flags for the 3 checkboxes

Drift acceptability limit

Pixel Size (µm)

String indicating side(s) where monolayer crosses the image